

EDAN



SonoTrax Series Ultrasonic Pocket Dopplers' superior sensitivity, high-fidelity sound, recording option and optional rechargeable battery ensure the doppler family to meet the routine examination requirements of obstetricians and midwives.

Specifications

- Display: 45 x 25 mm LCD display
- Heart Rate Range: 50~210 bpm
- Resolution: 1 bpm
- Recording Time: 240 s
- Audio Output Power: 0.5 W
- Auto Power-off: 1 minute with no signal or operation
- Size: 34 (D) x 89 (W) x 141 (H) mm
- Weight: about 300 g (including the battery)



Carry Case



Rolling Stand

SonoTrax Ultrasonic Pocket Doppler

- Easy-to-use
- Compact design
- Support 1.5 V AA battery
- Built-in audio recorder
- Large backlit LCD display
- Accurate fetal heart detection with high-fidelity sound
- Support 2 MHz, 3 MHz, 4 MHz, 5 MHz, 8 MHz waterproof probes
- Automatic and manual counting of FHR
- Interchangeable waterproof probe
- Automatic power-off



SonoTrax Lite



SonoTrax Basic



SonoTrax Vascular



SonoTrax Pro



SonoTrax II Pro